

UV PLUS 46G



UV PLUS 46G

- FOLDABLE INTRACULAR LENS WITH AN ADDITIONAL HEV (BLUE LIGHT) FILTER FOR PROPHYLAXIS OF AMD
- PROVEN HAPTIC DESIGN FOR SULCUS/SCLERA OR CAPSULAR BAG FIXATION

TYPE: FOLDABLE SULCUS / CAPSULAR BAG IOL

FEATURES: UV-PLUS -HEV FILTER, OPTIONAL SCLERA FIXATION.

OPTIC MATERIAL: BIOCOMFOLD

TOTAL DIAMETER: 15.0 MM

OPTIC DIAMETER: 6.0 MM

OPTIC: BICONVEX

INCISION SIZE: 2.5 MM WITH INJECTOR

A-CONSTANT: 119.1


ACD (ANTERIOR CHAMBER DEPTH): 5.61 MM


DIOPTR RANGE: +8.5 UP TO +35.0 DPT

ADVANTAGES

- PROTECT MACULA FOR UV AND HEV (HIGH ENERGY BLUE LIGHT) INCIDENCE OF LIGHT
- VERY GOOD BIOCOMPATIBILITY
- SMALL INCISION
- UNFOLDS SLOWLY AND CONTROLLED IN THE CAPSULAR BAG
- LOW PCO RATE
- MORE RESISTANT AGAINST YAG-LASER THAN ANY OTHER IOL MATERIAL

UV PLUS 46G





MORCHER® IOL



Morcher GmbH Kapuzinerweg 12 D-70374 Stuttgart

Tel.: +49 (0)7 11/9 53 20-0 Fax: +49 (0)7 11/9 53 20-80 e-mail: info@morcher.com www.morcher.com

Version:2006-09-01 Tsa